



Washington WIC – Certifier Competency Training Worksheet

6. Hematology

Competency	Certifier is able to:	Training Requirements
<p>Assess iron status and share with participant.</p>	<ul style="list-style-type: none"> • Use correct techniques to perform a blood test for adults, infants, and children. • Use correct techniques to perform non-invasive hemoglobin test for adults and children over two years, using Masimo Pronto device. • Enter hemoglobin in the participant’s file. • Share test results with participant and what it means. • State when to assess iron status for infants, children, pregnant, breastfeeding, and postpartum participants. <p>Note: Staff must be a Medical Assistant – Registered prior to performing blood tests.</p>	<ol style="list-style-type: none"> 1. Complete the following in the Learning Center (LC) <ul style="list-style-type: none"> • DOH STATE WIC Hematology Curriculum <ul style="list-style-type: none"> ○ For issues with images during the post test, contact: wawictraining@doh.wa.gov • Complete the post-test with 80% or higher score 2. Review all information below: <p>WIC Manual - Volume 1</p> <ul style="list-style-type: none"> • Chapter 10-Hematology • Chapter 14-Nutrition Risk Criteria <ul style="list-style-type: none"> ○ Review the Low Hemoglobin/ Hematocrit risk 201 <p>Core WIC Training</p> <ul style="list-style-type: none"> • Staff Tool: When to Test Hemoglobin or Receive Blood Work <p>Participant Nutrition Education Handout</p> <ul style="list-style-type: none"> • Pump up the iron

Be sure to discuss any outstanding questions or issues with certifier. Document completion of this competency on Certifier Competency Training Documentation Form.



Washington WIC – Certifier Competency Training Worksheet

6. Hematology

This institution is an equal opportunity provider.
Washington WIC does not discriminate.

To request this document in another format, call 1-800-841-1410.
Deaf or hard of hearing customers, please call 711 (Washington Relay)
or email wic@doh.wa.gov.



DOH 961-1118 May 2023

